

FRONT VIEW

PLUMBING CHASE

48" OR 58"

FOR ALL BUILT IN INSTALLATIONS
CALL TO CONFIRM DIMENSIONS

HOT AND COLD FIXTURE
FAUCET AND HOSE

37"

RETURN WALL

22"

SWING-IN SWING-OUT RAMF.

STRAINER

5"

WASTE
SUPPLY

SIDE VIEW

7"

HOT COLD

4"

21"

23"

OUTSIDE

28"

FINISHED FLOOR LINE

SWING-IN SWING-OUT RAMF.

NOTE: ADJUST LEVELERS FOR DRAIN PITCH

PET LIFT
EQUIPMENT CORPORATION®

VETLIFT
MEDICAL®

6340 17TH ST., CIRCLE EAST • SARASOTA, FLORIDA 34243

PLUMBING SPECIFICATIONS
AQUA QUEST 48" AQ 48K
AQUA QUEST 58" AQ 58 U

CUSTOMER SERVICE
TOLL FREE
888-388-1159